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MicroPort[®] REVIEW

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MicroPort Scientific Corporation
Registered Address: PO Box 309, Uglund House
Grand Cayman, KY1-1104, Cayman Islands

MicroPort Scientific Corporation
Business Address (HK): Level 54 Hopewell Centre,
183 Queen's Road East, HK

MicroPort Orthopedics Inc.
Add: 5677 Airline Road, Arlington,
TN 38002 Memphis, USA
Tel: +1 866 8720211
Fax: +1 855 4462247
Website: www.ortho.microport.com

MicroPort CRM, Clamart, France
Add: 4 Avenue R é amur, 92140 Clamart, France
Tel: +33 01 46 01 33 33
Fax: +33 01 46 01 34 58
Website: www.crm.microport.com

MicroPort CRM, Saluggia, Italy
Add: Via Crescentino s.n.c.,
13040 Saluggia (VC), Italy
Tel: +39 0161 4871
Fax: +39 0161 487681

MicroPort CRM, Santo Domingo,
Dominican Republic
Add: Zona Franca Industrial de Las Americas,
1-3 Km 22, Autopista Las Americas, Santo
Domingo, Dominican Republic
Tel: +809 549 9800
Fax: +809 549 9988

Medical Product Innovation, Inc.
Add: 199 Technology Drive, Suite 105,
Irvine, CA 92618, USA
Tel: +1 949 9883515
Fax: +1 949 9883517

MicroPort Orthopedics Japan K.K.
Add: Akasaka Enokizaka Building 10F, 1-7-1
Akasaka Minato-ku, Tokyo 107-0052, Japan
Tel: +81 3 6758 7279
Fax: +81 3 6758 7298

MicroPort Orthopedics Canada
Add: 3715 Laird Rd, Unit 9, Ontario, Canada
Tel: +1905 607 2929

MicroPort Scientific Coöperatief U.A.
Add: Paasheuvelweg 25, 1105 BP
Amsterdam, The Netherlands
Tel: +31 0 20 545 01 00
Fax: +31 0 20 545 01 09

MicroPort Scientific GmbH
Add: Harkortstr.11-13, 40880 Ratingen,
Germany
Tel: +49 2102 102 620
Fax: +49 2102 102 6210

MicroPort Orthopedics UK
Add: Willow House, Park West, Sealand Road,
Chester, CH1 4RJ, United Kingdom
Tel: +44 1244 572 120
Fax: +44 1244 572 165

MicroPort Orthopedics Srl.
Add: Via Liguria 18, 20068 Peschiera Borromeo
Milan, Italy
Tel: +39 02 51 69 91
Fax: +39 02 51 65 0279

MicroPort Orthopedics Brazil
Add: Rua Abid Auada, 35, Sao Paulo, Brazil
Tel: +5511 406 37171

MicroPort Orthopedics
Add: CCY1 UltraPark 1, Edificio/B Oficina
Heredia, 40104-Costa Rica
Tel: +5062 2099200

MicroPort Scientific India Pvt. Ltd.
Add: 801, 215 Atrium II, Behind Courtyard
Marriott, Andheri Kurla Road, Andheri (East),
Mumbai 400059
Tel: +91 22 71912900

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Position	Tissue Heart Valves	Mechanical Heart Valves	Surgical Annuloplasty Devices	Medtronic
1st	Medtronic	St. Jude Medical	Edwards Lifesciences	Medtronic
2nd	St. Jude Medical	Medtronic	Medtronic	St. Jude Medical
3rd	Edwards Lifesciences	LivaNova	St. Jude Medical	Edwards Lifesciences

Source: Millennium Research Group, Inc.

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Date	Product Type	Event
September 2017	TAVI	Venus A TAVI launched in China, first transfemoral device
August 2017	TAVI	J-Valve launched TAVI product, first transapical device
April 2016	TAVI	Venus Medtech (Hangzhou) acquires Transcatheter Technologies, a German manufacturer of transcatheter heart valves, Launched in 2018
April 2016	All Heart Valve Devices	Abbott Laboratories announces the acquisition of St. Jude Medical for \$25 billion (Michelle Cortez, 2016)
March 2016	Mechanical Heart Valves	Lepu Medical Technology (Beijing) expands its ownership of Lepu Pharmaceuticals from 60% to 100% (Reuters, 2016)
January 2016	Mechanical Heart Valves	Cryolife acquired On-X Life Technologies, a manufacturer of mechanical heart valves
October 2015	All Heart Valve Devices	Cyberonics merges with Sorin Group. The merged company rebrands itself as LivaNova
August 2015	TAVI Devices	Medtronic acquires Twelve for \$458 million (Medtronic, press release, August 25, 2015)
July 2015	TAVI Devices	Abbott Laboratories acquires Tendyne Medical for \$250 million (Abbott Laboratories, press release, July 30, 2015)
July 2014	TAVI Devices	Venus Medtech (Hangzhou)'s A-Valve device is among the first batch of products granted fast-track approval under the CFDA's Special Approval Process for Innovative Medical Devices
May 2012	TMVI Devices	Shanghai Zhongshan Hospital completes China's first three cases of TMVR using Abbott's MitraClip system

Source: Millennium Research Group, Inc.

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Product Type	Brand	Company	Launch Date	Status
TAVI Devices	Taurus One	PEIJIA	2021	Undergoing clinical trial
TAVI Devices	Cardioflow	Microport	2019	Waiting CFDA license approval
TAVI Devices	J-Valve	Jecheng (Suzhou)	2017	Estimate 350 cases commercialized in 2018
TAVI Devices	A-Valve	Venus Medtech (Hangzhou)	2018	Estimate 800 cases commercialized in 2018
TAVI Devices	CoreValve	Medtronic	Expected 2020	Received FDA approval in January 2014 Currently undergoing clinical trials in China
TAVI Devices	SAPIEN XT, SAPIEN3	Edwards Lifesciences	Expected 2020	Received FDA approval in June 2014 Currently undergoing clinical trials in China
TAVI Devices	MitraClip	Abbott Laboratories	Expected 2020 - 2030	Received FDA approval in October 2013 Currently undergoing preclinical studies in China
TAVI Devices	CardiAQ- Edwards	Edwards Lifesciences	Expected 2025 - 2030	FDA approval granted in 2015 to continue an early feasibility study

Source: Millennium Research Group, Inc.

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MicroPort Review: How long have you been using the product from MicroPort Orthopedics? Under what kinds of situation will you recommend patients to use MPO's products? And as for the young generation, and patients who ask for a strong function of sports, does medial pivot knee have any advantage for the patients?

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Dr. Lagast Joost: It has already been more than 20 years since I first used the Medial-Pivot Knee. I was the first user of medial pivot knee in the European Continent. I started on April 1998, and it was the first case on the EU. I have seen population that put in total knee replacement; the medial pivot is advanced than traditional designs, which has a much higher activity level. Because the traditional design were not regarding patients satisfactory or pain relief, but as you see 20 years ago the pain relief was the most important, but now we see the more active patients today are asking for more stability, which we obtain in the medial pivot. So for me, the future will be more and more important towards the medial pivot, we do not only want pain relief. The medial pivot was the first to copy the normal knee as close as possible; the other designs didn't look into the kinematics alignment. When the traditional knee only reach a 80% percent patient's satisfactory, medial pivot knee didn't randomly try other ways but to look into how normal knee works, then they look into the kinematics and discover the normal doesn't work in the way they thought, but in a medial pivot way. They transmitted the idea into the knee replacement. At first, the criticism was that the normal knee might work in a medial pivot way, but after the knee replacement, it will no longer move like that way afterwards. But the articles have proven that even after you have a total knee replacement, you still perform in the medial pivot way. You can always find new studying; the confirmation of the medial pivot is the future.

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MicroPort Review: You have just mentioned in the speech that Medial-Pivot Knee will be a “game-changer”, so can you explain the reason to the ordinary audience in several simple sentences why would it be a game-changer? As for MPO, we will keep evolving so that we can realize the ambition, thus for the view point of product itself, on what aspects should medial-pivot knee to improve further?

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Dr. Lagast Joost : It was really the game changer 20 years ago, when I first used it that many colleges said that it was crazy and told me not to believe it. But now you see that all companies are going towards all coming to the unequal lateral and external balancing, although they didn't look like the medial pivot, but they all coming towards the unequal medial and lateral components. Just proves that medial pivot works,

although they didn't admit 20 years ago, but they start to believe it is the right thing to do. But they are already too late, that's why I am glad to say I jump to the track first. Because I got the feeling it was the right thing to do and based on my philosophy based on my surgeries on the medial side and that why I was convinced at the first moment and now 20 years later, it is still the right answer. The first thing is not about innovation, it is to try to get it globally. To see how many surgeons are still not convinced by medial pivot, it is company's responsibility to show that everyone is going to the medial-pivot, but MicroPort is the only one that has 20 years clinical trial and hot high patients satisfactory we got. Let's see other companies they bring something new to the market that look like medial pivot, they don't have any follow up, and they still have to prove their idea. No company has brought a medial pivot product and proves it works and has brought that concept globally. Wright used to hold that part, but now it is for MPO you should push it more harshly and make it global that you have results that show it works, and you are still the pioneer. The revision system that comes has the same idea and again you can design the revision system that build on the same idea, as a pioneer.

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MicroPort Review:Despite the product itself, what should MPO improve on the marketing strategies so that MPO could expand its business landscape abroad? Compared with our competitors, what are our pros and cons?

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Dr. Lagast Joost: Try all big meetings, and introduce and promote the products with what you got like follow-up results. And I feel that still many surgeons they still don't know the medial pivot knee after so long time. If they start to look into the medial, it is up to company to convince them, in the long term, not short term. Other companies always go towards to the unequal medial and lateral compartment; because they all see normal knees they have more stability in the medial side than the lateral side, and they try to introduce those things, and the philosophy behind is that you need to improve the stability on the medial side, the symmetrical design is not really like it should be, thus they should update the old design. For many years, I have been done mechanics, compared with kinematic, you need to do less releases for kinematics, for sure you do mechanics. For the medial pivot, you don't have to equal balancing, you go to unequal balancing, you have the more stability on the medial side and looser on the lateral side, that the key issue for the kinematics alignment, if you go unequal, it is the key issue. If you need to do the normal knee, you need to copy the normal knee. I always try to follow this rule.

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MicroPort Review:The surgical observation visitation in Belgium hold in June left us a deep impression. A group of domestic and abroad surgeons benefited a lot from it. How do you evaluate the value of such visitations? On what aspects could MPO improve and update in the future?

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Dr. Lagast Joost : First of all, you have to adapt a little bit on the unequal balancing. You must have the surgeons convinced by your idea, so I am convinced and I want to start. At the beginning of the study, you have a good selection of surgeons. Absolutely it is a must they go for training before you start, or else they would do it in the old way. You have to do it, and you will start the medial pivot ever since. I have seen results could be not good, but it's not due to product, but due to the training. To see the surgeons coming to training, since they got convinced. There are always surgeons told me that they felt lucky to be there, so that they can learn the skills in that way, and it really works. They should start the surgeries with fully training, and understand the concept of the medial pivot.

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MicroPort Review:What do you think of this year's COA? Is there any advanced theory or tech that gets your attention? What will be the direction for the tech to innovate in the orthopedics area? What kinds of opportunities and challenges will these innovations mean to MPO in the coming future?

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Dr. Lagast Joost : For the moments there are something that absolutely catches me, the major point for the discussion should be robot or not robot. Be honestly, I am still not convinced, because if you look in to the eye of navigation, you could see how many surgeons still using the navigation. You will see 90% of the surgeons no longer using it, and I think it will be a little bit same for the robot. I do say that you still need to open to see the new tech including in the surgery, but the main concept will still be medial pivot. The product will never be perfect, but there will be new tools for us to improve the result of patients. "Patients always come first" is what we have to keep in mind. If any new tech can improve the result of the surgeries, we should look at it. New material, for example the polymer, it will always be an issue in the long term. If there are some new materials that will stand for 40 years, we will no longer need revisions. That would be the idea.

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MicroPort Review: Have you been to China and did communications with Chinese orthopedic surgeons? What would you think of the Chinese surgeons, their techniques and their philosophy of joint replacement?

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Warwick J.M. Bruce: Yes, by the operating and lectures. Every time I came to China, there's an exponential improvement in the knowledge of the surgeons. The standard of the hospitals has improved every single time. I found I have never been disappointed in the skills of Chinese surgeons. I've been to the most western city in China and the skill of the surgeons there was amazing. So I think the only thing that was holding back the Chinese surgeons in the past was that they didn't have all the options available to do difficult joint replacements, including instruments and prosthesis. But everything that I was noticing every time, it got much better.

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The other thing that's improved, which is helpful, is that the Chinese surgeons' English have improved so much that they have a wide education sphere. So they can look up European literature that's written in English or American literature, whereas in the past they mostly did in Chinese, especially for young generations. Just like I will read European articles in English and American articles, I get different philosophies.

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MicroPort Review: Do you think there is any gap between Chinese surgeons and Australian surgeons? What's that?

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Warwick J.M. Bruce: I think in the major cities there is no gap, while in the lesser cities; the second tier and third tier, there are some gaps. But it's not just in surgical technique education. It's because there's equipment, education, things like these. And I think if you go to any country that has a huge population and second and third tier cities, all the best-educated and most experienced people are in the first tier. So I think it's just the China has to embrace education from first tier cities, educating second and third tier cities, and for international education. Just like we have international education, Americans, English, French come to Australia to talk. You will get different philosophies from different countries. It's the responsibility that the surgeon from tier one cities, they go to other cities and express the education. If I was running education in China, I'd have first tier of orthopedic surgeon cities but with other international, so they would come and do meetings with Chinese surgeons from first tier institutions and international surgeons, like a user group.

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MicroPort Review: In your knowledge, what a good knee prosthesis design features? What kinds of facts do they have to promote to go further to get the good outcome of the patients?

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Warwick J.M. Bruce: Firstly, it's the surgical technique that the most important in all designs. Secondly, number one knees from the number one companies do well, but

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there are design parameters that improve, particularly knee replacement. In hip replacement, it's less obvious. So certainly, the medial pivot knee has theoretical advantages. It has experimental advantage like fluoroscopy and it has clinical advantages. And I think that is a knee that has an advantage, but not if it's not done properly. It is all about education and experience. So, certainly the medial pivot knee for one requires slightly different philosophy on how to do on the knee, as we've talked in the past. So it's very important that surgeons visit people, do the particular case to see it, and do category workshops if possible. Assist surgeons who do it, and then do their own.

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MicroPort Review: What's your expectation for the patient outcome with the Medial Pivot Knee? For the patients' outcome, the relief of pain, ranges of motion and something else?

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Warwick J.M. Bruce: I think patient outcome always improved with education and good design, and in the medial pivot principle, obviously a very good design. If you do good surgeon education, you get a very good result. Regarding the relief of pain, I think this is where western and eastern differing. So in western patients, they are very demanding of function. And the reason they demanding of function is because they don't have that much pain before the operation.

Because of obesity, we are replacing knees that don't have a lot of arthritis, but because they take a lot of blood. So it's very difficult in western society because patients are out of condition and overweight. In Asian societies, the patients are not overweight, but they have significant deformity. And they have a lot of pain. So I think in Asian society, it's about relieving pain and function is second. However, functions very important for younger patients and for people who have higher expectations.

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MicroPort Review: The Medial Pivot Knee has been twenty years clinical trial. So compared to the conventional prosthesis, what was the advantage?

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Warwich J.M. Bruce: The advantage is what the patients feel. Firstly, we have no problem that lasts as long as other knees. We have no problems that it wears better than other knees or equally. We have no problem that bends as well about other knees. But we certainly notice that functionally the patients have less instability. They're better at stair walking. They're better at kneeling and they better at crouching.

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MicroPort Review: How long you performed on the medial pivot knee and what's the outcome of patience? What's your experience on medial pivot knees?

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Warwich J.M. Bruce: I did a series of medial pivot knees nineteen years ago and there have been no failures and one reoperation for a patella button, so there is no doubt it's in the top echelon of results. The issue I noticed with the patients both in the medial pivot knee and other designs of knees that incorporate a medial pivot principle, they are better with stairs. They have less AP drops and they have good function.

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MicroPort Review: MicroPort Orthopedics is devoted to events or develops some Chinese or Asian customized prosthesis. What aspects do you think for us manufacturer to promote?

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Warwich J.M. Bruce: Firstly, I think MicroPort should be very heavily invested in education. Secondly, because you have purchase the medial pivot design from an American company, I think you should look seriously at sizing. So maybe you have one or two smaller sizes and maybe a little

bit thinner for mediolateral. But that would just require study looking at CT scans of Asian knees in China and maybe tweaking the design a little bit smaller sizes may be a little bit thinner from medial to lateral. When I did studies in the past, for western knees that being applied to eastern knees, I'd like to have them to be smaller and thinner.

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MicroPort Review: How do you think of the robot being applied to the TKA and THA? Can you predict the trends for the joint replacement?

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Warwich J.M. Bruce: The advantage of the robot is that you have reproducible information. When you start the knee, it opens one millimeter or one degree medial two degrees lateral. You know exactly when you finish the operation, it open one medial two lateral. And then you can retrospectively look how the knee performs over time, and know maybe that's too loose or too tight. The other advantage of the robot is you have combined gap balancing and mechanical alignment and kinematics, so you can combine both very accurately. However, I think the robot is less important than spending money on education in China, because you have a lot of people that made knee replacement. If you train surgeons well, they're going to do a good job. So that's way. I think you should choose the best bang for your buck is. However, the advantage of the robot is that it is a very good marketing tool. That's more important to companies, less important to me. I think there are a couple of trends that I can predict. I think that robots

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will become more important, especially for knees. I think preoperative analysis of pelvic rotation is very important for hips. I think also the companies will play a greater role in the clinical delivery of the knee replacement. And that means to me, and this is what I've always thought should happen in western society, is that the company supplies the assistor, who is also the company representative, so that there's less chance of things going wrong. But of course in Australia, the company rep is the teachers, they assist everything so if you're going to do it that way, you might as well have the company rep doing the operation with the surgeon. That's one trend. I think the other trend is that we will make the polyethylene so good that they won't wear out.

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MicroPort Review:2018 is the twentieth year anniversary for medial pivot knee and for MicroPort. Do you have any wishes to us? As a company like a MicroPort, we are devoted to promoting a professional education through internet that connects everyone; Do you have any advice on personal education?

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Warwick J.M. Bruce:Of course! I wish you the best I think the other important thing is that MicroPort is a Chinese

company. It's very well respected in China, which is very important. So they have American manufacturing and Chinese know how. And that's very important for an international reputation, because western society has still probably not come to realize the sophistication of the Chinese manufacturing. They trust American manufacturers because that that's what they've had. But they also know Chinese innovation and the combination is extremely important for an international market. You only have to look at what MicroPort has done in cardiac to know it's a force to be reckoned with in the world. And it got good products. It got innovation from China. It got all will charm from America and manufacturing. So it's going to be a very successful company.

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
Tubridge[®]
血管重建装置
Vascular Reconstruction Device

Parent Artery Reconstruction for Large or Giant Cerebral Aneurysms Using the Tubridge Flow Diverter: A Multicenter, Randomized, Controlled Clinical Trial (PARAT)

X J.-m. Liu, X Y. Zhou, X Y. Li, X T. Li, X B. Leng, X P. Zhang, X G. Liang, X Q. Huang, X P.-f. Yang, X H. Shi, X J. Zhang, X J. Wan, X W. He, X C. Liang, X G. Zhu, X Y. Xu, X B. Hong, X X. Yang, X W. Bai, X Y. Tian, X H. Zhang, X Z. Li, X Q. Li, X R. Zhao, X Y. Fang, and X K. Zhao, for the PARAT investigators

American Journal of Neuroradiology May 2018, 39 (5) 807-816

ABSTRACT

BACKGROUND AND PURPOSE: Although flow diverters have been reported with favorable clinical and angiographic outcomes in various literatures, randomized trials determining their true effectiveness and safety are still in lack. The Parent Artery Reconstruction for Large or Giant Cerebral Aneurysms Using the Tubridge Flow Diverter (PARAT) trial was designed to evaluate the safety and efficacy of the Tubridge flow diverter in the treatment of large or giant aneurysms in comparison with Enterprise stent-assisted coiling.

MATERIALS AND METHODS:

This prospective, multicenter, randomized trial was conducted at 12 hospitals throughout China. Enrolled adults with unruptured large/giant intracranial aneurysms were randomly assigned (1:1) to receive either Enterprise stent-assisted coiling or Tubridge flow diverter implantation. The primary end point was complete occlusion at 6-month follow-up, while secondary end points included technical success, mortality, target vessel-related stroke, aneurysm bleeding, in-stent stenosis, parent artery occlusion, and the frequency of all adverse events.

RESULTS:

Among 185 enrolled subjects, 41 withdrew before procedure initiation. Overall, 82 subjects underwent Tubridge implantation, and 62 subjects were primarily treated with stent assisted coiling. The results of 6-month follow-up imaging included complete occlusion rates of 75.34% versus 24.53% for the Tubridge and stent-assisted coiling groups, respectively, with a calculated common odds ratio of 9.4 (95% confidence interval, 4.14 – 21.38; $p < .001$). There was a higher, nonsignificant frequency of complications for Tubridge subjects. Multivariate analysis showed a decreased stroke rate at the primary investigational site, with a marginal P value ($p = .051$).

CONCLUSIONS:

This trial showed an obviously higher rate of large and giant aneurysm obliteration with the Tubridge FD over Enterprise stent-assisted coiling. However, this higher obliteration rate came at the cost of a nonsignificantly higher rate of complications. Investigational site comparisons suggested that a learning curve for flow-diverter implantation should be recognized and factored into trial designs.

ABBREVIATIONS:

FD: flow diverter; LTF: lost-to-follow-up; PARAT: Parent Artery Reconstruction for Large or Giant Cerebral Aneurysms Using the Tubridge Flow Diverter; SAC: stent-assisted coiling.

Received January 16, 2018; accepted after revision February 26.

From the Department of Neurosurgery, (J.-m.L., Y.Z., P.-f.Y., B.H., Y.X., Q.L., R.Z., Y.F., K.Z., Q.H.), Shanghai Changhai Hospital, Shanghai, China; Department of Neurosurgery (Y.L., X.Y.), Tianan Hospital, Beijing, China; Department of Interventional Radiology (T.L., W.B.), Henan Provincial People's Hospital, Henan Province, China; Department of Neurosurgery (B.L., Y.T.), Huashan Hospital, Shanghai, China; Department of Neurosurgery (P.Z., H.Z.), Xuanwu Hospital, Beijing, China; Department of Neurosurgery (G.L., Z.L.), General Hospital of Shenyang Military, Shenyang Province, China; Department of Neurosurgery (H.S.), First Clinical Hospital affiliated with Harbin Medical University, Heilongjiang Province, China; Department of Neurosurgery (J.Z.), Second Affiliated Hospital of Zhejiang University, Hangzhou Municipality, Zhejiang, China; Department of Neurosurgery (J.W.), Renji Hospital, Shanghai, China; Department of Neurosurgery (W.H.), Second Affiliated Hospital of Guangzhou Medical University, Guangzhou, China; Department of Neurosurgery (C.L.), First Hospital of China Medical University, Shenyang

Municipality, Liaoning, China; and Department of Neurosurgery (G.Z.), Southwest Hospital, Chongqing Municipality, China.

Jin-min Liu and Yu Zhou contributed equally to this work.

The PARAT investigators are listed in the On-line Appendix.

The trial was funded by the National Research and Development Project of Key Chronic Diseases (Grant No. 2016YFC1300700), the National Science and Technology Supporting Program (2011BAI08B14), and the project of the Shanghai Science and Technology Commission (11DZ1921603) and was sponsored by unrestricted grants from MicroPort NeuroTech Co, Shanghai, China.

Please address correspondence to Jian-min Liu, Department of Neurosurgery, Changhai Hospital, Second Military Medical University, 168 Changhai Rd, Shanghai 200433, China; e-mail: liu118@vip.163.com

Indicates open access to non-subscribers at www.ajnr.org

Indicates article with supplemental on-line appendix and tables.

Evidence-Based Medicine Level 1.

<http://dx.doi.org/10.3174/ajnr.A5619>

Large (10 – 25 mm) or giant (≥ 25 mm) aneurysms are a very challenging subtype among intracranial aneurysms, with much higher risks of rupture and poorer outcomes compared with small aneurysms.^{1,2} Despite recent technical advances, the treatment of large or giant aneurysms remains technically disappointing, with high complication and/or recanalization rates.³⁻⁶

In recent years, various flow diverter (FD) devices, such as the Pipeline Embolization Device (Covidien, Irvine, California), the Silk flow diverter (Balt Extrusion,

Montmorency, France), the Flow-Redirection Endoluminal Device (FRED; MicroVention, Tustin, California), and the Surpass stent (Stryker Neurovascular, Kalamazoo, Michigan), have been increasingly used in >50 countries. These devices were believed to improve long-term effectiveness, due to their capability to alter intrasaccular hemodynamics.⁷ Unfortunately, reported clinical results have varied significantly, with aneurysm occlusion rates and periprocedural complication rates ranging from 49% to 93.4% and 2.8% to 11%, respectively. On the other hand, there is only 1 randomized controlled trial reported until now (Flow Diversion in Intracranial Aneurysm

Treatment [FIAT]), which showed very high complication rates and below-expectation effectiveness.⁸ These results confused neurointerventionalists about the true effectiveness and safety of FDs. The recent publication of Raymond et al⁸ echoed these sentiments by suggesting that more randomized trials are needed to determine the role of flow diversion in the management of aneurysms.

In 2012, a lack of prospective data and reports of significant adverse outcomes associated with FD use in treating intracranial aneurysms fueled our interest in leading a multicenter, prospective, randomized, controlled trial (Parent Artery Reconstruction for Large or Giant Cerebral Aneurysms Using the Tubridge Flow Diverter [PARAT]) assessing FD treatment of large and giant intracranial aneurysms. In contrast to the FIAT trial, the PARAT trial focused on unruptured large or giant internal carotid artery or vertebral artery aneurysms, the traditional indication for FD implantation. By comparing it with a well-established treatment, stent-assisted coiling (SAC), we attempted to characterize the safety and effectiveness of the Tubridge FD (MicroPort Neuro-Tech, Shanghai, China) in this specific subset of intracranial aneurysms.

MATERIALS AND METHODS

Study Design and Participants

The PARAT trial used a prospective, multicenter, parallel-group design, with balanced randomization (1:1). The primary trial purpose was to compare outcomes of subjects with unruptured large/giant intracranial aneurysms who were treated with either Enterprise SAC (Codman & Shurtleff, Raynham, Massachusetts) or Tubridge FD implantation. The study was conducted in 12 centers throughout China, each of which was required to have performed >30 stent-assisted coiling procedures per year during the previous 5 years. The PARAT trial protocol was approved by all relevant local ethics boards.⁹ Site investigators generated trial data, with monitoring and data base maintenance completed by a commercial clinical research organization. The corresponding author had full access to all trial data and had final authority for key decisions relevant to this publication. At each trial center, prospective subjects were screened for trial eligibility, based on having an unruptured ICA or vertebral artery saccular aneurysm (including recanalized

aneurysms) measuring ≥ 10 mm in maximum diameter and ≥ 4 mm across the aneurysm neck. Investigators recruited only those who met all inclusion criteria and none with the exclusion criteria. Patients with ruptured aneurysms or other intracranial diseases were excluded. Detailed inclusion and exclusion criteria are listed in On-line Table 1. Investigational sites did not record or retain a log of patients who were screened for eligibility. Written informed consent was obtained from each participant before enrollment.

Randomization and Masking

After we obtained each prospective subject's informed consent, treatment allocation was initiated by a researcher contacting a clinical research associate who was independent of the patient recruitment procedure. Subject randomization was accomplished through an interactive Web response system, which was developed by information technology managers from an independent clinical research organization under the instruction of a clinical trial statistician. The randomization scheme included built-in stratification by the participating center, aneurysm size (≤ 15 mm versus >15 mm), and aneurysm location (anterior circulation versus posterior circulation). The on-line central randomization ensured that the allocation sequence was concealed from investigators who were recruiting patients before the decision to randomize. Treatment groups were balanced using minimization criteria. Masking of the randomization results to local investigators or trial participants was impossible to achieve. However, members of an independent Core Laboratory and those who followed up with trial subjects via telephone or face-to-face on-site interviews were blinded as to treatment allocation.

Interventions, Treatment, and Evaluation Procedures

For all recruited patients, dual-antiplatelet drugs (300 mg/day of aspirin plus 75 mg/day of clopidogrel) were prescribed for at least 3 days before the pivotal procedure. All FD placement procedures were performed with the patient under general anesthesia and via a transfemoral approach. After sheath placement, heparin (100 u/kg) was administered to maintain an activated clotting time of 250 – 350 seconds throughout the procedure. Next, a suitable guiding catheter (7F for subjects

Table 1: Structural comparison between different FDs

FD Type	Size	Braided Microfilaments ^a		Radiopaque Microfilaments	Flared End	Metal Coverage	Retrievable
		No.	Material				
Pipeline	3–5.5 mm	48	75% cobalt chromium and 25% platinum	NA	No	30%–35%	Yes ^b
Silk	2.5–5 mm	48	Nickel-titanium alloy	4 Platinum wires	Yes	35%–55%	Yes
Surpass	2.5–5 mm	2.5 mm (36) 3 and 4 mm (60) 5 mm (84)	Cobalt-chromium	12 Platinum wires	No	30%	NA
FRED	2.5–5.5 mm	Inner layer (48) Outer layer (16)	Nickel-titanium alloy	2 Tantalum wires	Yes	NA	Yes
Tubridge	2.5–6.5 mm	3.5 mm (46) 3.5 mm (62)	Nickel-titanium alloy	2 Platinum-iridium wires	Yes	30%–35%	Yes

Note:—NA indicates not applicable.

^a Braided microfilaments in this table mean those main wires excluding microfilaments especially for radiopaque usage.

^b Pipeline Flex embolization devices are retrievable.

to be treated with Tubridge and coils) was placed in the distal internal carotid or vertebral artery. Subsequently, treatment group subjects underwent Tubridge FD implantation alone or in combination with bare coils, while control group subjects were treated by Enterprise stent-assisted coiling with bare coils. The Tubridge Flow Diverter, developed by MicroPort Neuro-Tech, Shanghai, China, is designed to encourage the formation of intra-aneurysmal clot, while concurrently repairing and reconstructing the parent artery. Previously, we have characterized the Tubridge FD and its structural differences from the Pipeline and Silk flow diverters, such as the platinum-iridium material used for the radiopaque microfilaments, more braided microfilaments for the large-sized FD, and a decreased incidence of shortening (Table 1).⁹ Deployment of the Tubridge FD requires a combination “push” and “pull” technique to ensure full opening of the device and to increase the metal coverage across the aneurysm neck, as previously described.⁹ The approach and timing of coil placement and the number of Tubridge devices implanted were left to the discretion of each operator. For control subject aneurysms, Enterprise stents were implanted according to the instructions for use of the product, and all aneurysms were embolized (coiled) to achieve maximum packing density. All subject treatments were well-documented, and details included aneurysm shape and dimensions, width of the aneurysmal neck, diameter of the parent artery, and all relevant treatment outcomes. The postoperative antiplatelet regimen was as follows: <6 weeks: 300 mg of aspirin +75 mg of clopidogrel; 6 weeks to 3 months: 100 mg of aspirin +75 mg of clopidogrel; and, ≥ 3 months: 100 mg of aspirin indefinitely.

Digital copies of angiograms, including 3D rotational angiography as well as working projection images and other anteroposterior/lateral angiographic images, were collected by the clinical research organization – assigned clinical research associate of each site. These were sent to an independent Core Laboratory for analyses by up to 3 experienced neurointerventionalists (“analysts”). On the first pass, 2 analysts separately reviewed each subject's imaging. If the initial 2 interpretations conflicted, a third analyst provided a tie-breaking assessment. Follow-up angiographic results were classified into 4 categories, according to the immediate degree of embolization: 1) occluded, defined as no contrast filling into the aneurysm sac; 2) improved, defined as decreased contrast filling into the aneurysm sac; 3) stable, defined as unchanged contrast filling into the aneurysm sac; and, 4) recanalized, defined as increased contrast filling into the aneurysm sac. Analysts also documented instances of parent artery occlusion or in-stent stenosis.

At 1, 3, and 6 months postimplantation, experienced site investigators followed up with each subject by telephone or by in-person clinical interviews. To provide some level of objectivity, these interviewers were blinded to the treatment allocation. If any adverse events were identified, investigators documented, in full detail, the following information: symptoms, event duration and severity, possible causes and associations, actions taken, and event resolution/final outcomes. A fully independent Clinical Events Committee evaluated all investigator-documented adverse events and categorized them according to the most likely causal relationship: procedure-

related, device-related, disease-related, or unrelated to either disease, device, or procedure.

Study Outcomes

The primary end point was complete aneurysm occlusion at the 6-month follow-up. During reviews of 3D rotational angiography and other angiographic images, members of the Core Laboratory consistently followed precise analytic criteria for establishing the final effectiveness end point. Only those aneurysms exhibiting no contrast filling in the aneurysm were judged as meeting the complete occlusion designation. End point analysis was expressed as percentage occurrence within the test and control groups.

Secondary end points included the following:

- 1) Immediate technical success rate, which included successful device delivery, exact stent positioning, and full expansion of the devices.
- 2) Death or stroke related to target vessel (assessment time points: 30 days, 90 days, and 1 year postoperation). Stroke was defined as sudden symptoms and signs of a focal cerebral function deficit associated with cerebral circulation disorders and including hemorrhagic or ischemic stroke. Hemorrhagic stroke was defined as acute extravasation of blood into the brain parenchyma or subarachnoid space with associated neurologic symptoms, whereas ischemic stroke was defined as rapid onset of a new focal neurologic deficit or rapid worsening of an existing focal neurologic deficit with clinical evidence of infarction not attributable to a nonischemic etiology (not associated with brain infection, trauma, tumor, seizure, severe metabolic disease, or degenerative neurologic disease).
- 3) Aneurysm bleeding rate (assessment time points: 30 days, 90 days, and 1 year postoperation), including intraoperative rupture and delayed aneurysm rupture confirmed by CT.
- 4) The rate of in-stent stenosis (assessment time point: 6 months postoperation).
- 5) The rate of parent artery occlusion (assessment time point: 6 months postoperation).
- 6) General adverse events (assessment time points: 30 days, 90

days, and 1 year postoperation).

Statistical Analysis

On the basis of an assumption of a 20% lost-to-follow-up (LTF) rate, a 124-subject sample size was planned (62 per group) to obtain a valid result (2-tailed test, significance level of $\alpha=.05$, power of $1-\beta$ 0.80). However, because postrandomization subject withdrawals were higher than anticipated, the independent Data Monitoring Committee of the trial suggested a larger sample size to reach the required minimum number for each group (62 subjects). Differences between the PARAT trial protocol and a final report were included in the statistical analysis plan.

Data analyses included a modified intention-to-treat approach based on a population of subjects who were recruited and treated (full analysis set). Those who left the trial immediately after randomization and who did not undergo the pivotal procedure were removed from this analysis. Those who completed treatment and the entire follow-up protocol constituted the perprotocol set, while those who were treated and had at least 1 safety evaluation were included in the safety set. To verify intergroup balance, the statisticians compared the baseline characteristics of subjects between the 2 study arms. Proportions were used for categorical variables, and medians with interquartile ranges were used for continuous variables. The categorical variables were compared between study arms using the χ^2 or Fisher exact test. For continuous variables, the t test or Mann-Whitney rank test, as applicable, was used according to the distribution of the data.

Analysis of the primary effectiveness end point involved a comparison of complete occlusion rates at the 6-month follow-up between the treatment and control groups in the full analysis set using the χ^2 method in combination with multivariable logistic regression adjusted for subject age (60 years or younger versus older than 60 years) and aneurysm size (≤ 15 versus >15 mm). A sensitivity analysis was performed to assess the effect of missing data (ie, those who were partially or wholly LTF after undergoing the pivotal procedure). A per-protocol set analysis was also performed excluding the following: subjects who were enrolled in the trial but who did not undergo the index procedure, those in whom the procedure failed, and those who

did not reach the primary end point because they were deemed LTF. Secondary outcomes were compared between groups in the safety set using the χ^2 or Fisher exact test, as applicable. An additional multivariable logistic regression was performed to explore factors potentially affecting subject outcomes.

Statistical analyses were conducted by statisticians at the Institute of Clinical Evaluation, affiliated with Beijing University, and the data were analyzed using Statistical Analysis System software, Version 9.2 (SAS Institute, Cary, North Carolina). All tests were 2-sided, and a p value $< .05$ was considered statistically significant.

An independent Data and Safety Monitoring Committee had unrestricted access to trial data, to allow periodic monitoring of trial progress. This trial was registered on the Chinese Clinical Trial Registry: ChiCTR-TRC-13003127.

RESULTS

Subject Characteristics and Disposition

Between December 18, 2012, and May 9, 2014, one hundred eighty-five subjects were enrolled and randomized. The distribution of subjects among the investigational sites is shown in On-line Table 2. Of considerable surprise and prompting a midstudy revision of enrollment strategy was the high proportion of subjects ($n=41$) exiting the study after randomization. Details of these LTF subjects are found in the final part of this article in a discussion of trial limitations as well as in the Figure.

Pivotal treatments were initiated in the remaining 144 subjects who formed the full analysis set based on a principle of modified intention-to-treat. Due to a tortuous parent artery or a wide aneurysm neck, investigators failed to catheterize the parent artery or deliver devices across the aneurysm neck in 6 instances (3 per trial group). Of 138 treated subjects, 5 died; 1 withdrew after serious procedure-related complications; and 6 were defined as LTF without reason. The remaining 126 subjects were included in the per-protocol analysis. There were no subject crossovers among trial groups. Treatment and follow-up details are shown in On-line Table 3.

A comparison of baseline characteristics between the 2 groups in the full analysis set showed similar distributions except for subject

age ($p=.036$, Table 2), which was not considered a prognostic factor for the primary outcome. Because a high proportion of participants quit after randomization, we also compared baseline characteristics between the 41 excluded subjects and the 144 subjects in the full analysis set and determined that the distributions were comparable (On-line Table 4).

Primary Outcome

As shown in Table 3, in the 6-month angiographic follow-up, aneurysms treated with Tubridge FDs were associated with a favorable shift toward a complete occlusion rate. In the per-protocol analysis of 126 cases for the primary effectiveness end point, the complete occlusion rate was 75.34% (55/73) for the Tubridge group and 24.53% (13/53) for the Enterprise control group, with an adjusted common odds ratio of 9.31 (95% confidence interval, 4.00 – 21.66; $p<.001$). Unadjusted results and sensitivity analyses also showed superior results for the Tubridge group. Table 4 lists 6-month 100% aneurysm occlusion rates in the context of aneurysm location and size.

Secondary Outcomes

Immediate technical success rates did not differ significantly between the 2 trial groups. Three procedures in each group failed due to difficulties in advancing a microcatheter into the distal arteries. The calculated technical success rates were 96.34% (79/82) and 95.16% (59/62) in the Tubridge FD group and Enterprise control group, respectively, with a common odds ratio of 1.34 (95% CI, 0.26 – 6.87; $p=.726$).

Trial statisticians calculated a trend toward an increased but nonsignificant risk of complications for subjects in the Tubridge group. Hemorrhagic stroke occurred in 7 subjects, including 5/82 (6.1%) in the Tubridge group and 2/62 (3.23%) in the control group. Ischemic stroke, related to target vessels, occurred in 8/82 (9.76%) of the Tubridge group and 6/62 (9.68%) of the control group. These complications resulted in 4 deaths: Three were Tubridge subjects and 1 was an Enterprise-treated subject. In addition, 1 subject in each of the 2 trial groups died during follow-up for reasons unrelated to target vessels. Overall, the rates of death or stroke related to target vessels in 1-year follow up were 14.52% and 17.07% in the control and Tubridge groups, respectively. Statistically, there was no significant difference between the groups ($p=.678$).


FIGURE. Flow diagram of the PARAT trial.

The aneurysm bleeding rate was comparable between groups, with 3.66% versus 1.61% for the Tubridge group and Enterprise control group, respectively ($P=0.634$). Overall, general adverse events occurred in 46/82 (56.10%) of the Tubridge group and 33/62 (53.23%) of the Enterprise control group; among these, most were symptoms unrelated to the device, index procedure, or aneurysmal disease, such as headache, vomiting, or fever and there was no significant difference between the groups at the 30-day, 90-day, or 1-year follow-up. There was a trend toward a higher rate of in-stent stenosis or parent artery occlusion in the Tubridge group, but the difference was not statistically significant (Table 5).

Additional Complication Details

As shown in On-line Table 5, target vessel – related stroke, both ischemic and hemorrhagic, occurred in 21 subjects. Eight occurred during the procedure; 9, within the first 30 days after

the index procedure (procedure-related); and the remaining 4 cases, at 1, 2, 3, and 7 months posttreatment.

Intraoperative bleeding occurred in 3 subjects. One Tubridge subject had a microwire injury to a distal artery and died, while the other 2 subjects with intraoperative bleeding (1 per group/mRS of 1 at follow-up) recovered well after prompt coil deployment immediately after the onset of rupture. The Clinical Events Committee concluded that these complications, though unrelated to the Tubridge and Enterprise devices per se, were classified as procedure-related complications.

Periprocedural bleeding occurred in 4 cases (3 Tubridge/1 Enterprise). Delayed aneurysm rupture accounted for 2 complications (both large or giant aneurysms treated with FD alone). Of the remaining 2 cases, 1 Enterprise subject presented with cerebral hematoma and 1 Tubridge subject presented with

mild subarachnoid hemorrhage in the Sylvian fissure. Both were thought to be unrelated to aneurysm rupture. The patient with a cerebral hematoma died due to prehospital delay, while the other one recovered well.

Overall, there were 14 cases of target vessel – related ischemic stroke (8 Tubridge/6 Enterprise). Five subjects experienced intraoperative ischemia; another 5, within the first 30 days; and the remaining 4, at 1 – 7 months posttreatment.

Among those with intraoperative ischemic events, 3 were Tubridge subjects. One had a stroke after anterior choroidal artery occlusion due to coil introduction; another had a frontoparietal acute infarction due to parent artery occlusion caused by intrastent thrombosis; and, the other one was thought to be the result of a thromboembolic event. All of them resulted in contralateral limb weakness, and symptoms resolved in 1 case. In the 2 control subjects, patients had ischemic symptoms after uneventful procedures. Thromboembolic events were considered after further CT and angiographies, and both recovered well after medical therapy.

Among 5 subjects with periprocedural thromboembolic events (2 Tubridge/3 Enterprise), ischemic symptoms were all identified within 1 week after the procedure. Symptoms resulting from these periprocedural ischemia events were characterized as mild and transient, with all events resolving well after aggressive antiplatelet therapy.

There were no hemorrhagic events after the periprocedural period. However, delayed ischemic stroke occurred in 4, including 3 Tubridge subjects and 1 Enterprise subject. For the 3 Tubridge subjects, 2 strokes were identified with parent artery occlusion after further examination and the patients were treated conservatively; and the other 2 patients (1 per group) had ischemic symptoms 1 and 2 months postprocedure, with patent parent arteries. All these symptoms improved during the follow-up period.

Overall, among all these complications, there were 4 deaths and 2 disabilities identified during the last follow-up. The procedure-related mortality rate was 3.66% (3/82) in the Tubridge group versus 1.61% (1/62) in the Enterprise group, while the procedure-related morbidity rate was 2.4% (2/82) in the Tubridge group versus zero in the control group in 1-year

follow-up.

Multivariate Analysis

To further explore the factors affecting angiographic outcomes and complication occurrences, trial statisticians performed multivariate analysis, including subject age, aneurysm size, aneurysm location (anterior circulation versus posterior circulation), and treatment technique (FD alone; FD + coils; and Enterprise stent-assisted coiling). The primary investigational site (Shanghai Changhai hospital) showed a decreased stroke rate compared with other trial centers, with a marginal P value ($OR=4.81$; 95% CI, 0.991 – 23.335; $p=.051$).

DISCUSSION

Since the introduction of flow diverters to the global neurovascular community, many clinical studies were conducted to evaluate the safety and effectiveness of FDs, such as the Pipeline for Uncoilable or Failed Aneurysms (PUFS), Aneurysm Study of Pipeline in an Observational Registry (ASPIRe), and PITA trials.^{8,10-18} Additional effort to evaluate the effectiveness and safety of FDs versus conventional treatments included the Multicenter Randomized Trial on Selective Endovascular Aneurysm Occlusion with Coils versus Parent Vessel Reconstruction using the SILK Flow

Table 2: Baseline characteristics of the full analysis set

Characteristics	Tubridge Group n = 82	Enterprise Group n = 62
Mean age (yr)	52.11 (10.31)	55.66 (9.53)
Sex ratio (male/female)	21:61	13:49
Medical history		
Cerebrovascular stroke	9 (10.98%)	8 (12.90%)
Coronary artery disease	1 (1.22%)	4 (6.45%)
Hypertension	41 (50.0%)	25 (40.32%)
Hypercholesterolemia	4 (4.88%)	3 (4.84%)
Diabetes	5 (6.10%)	2 (3.23%)
Current/previous smoking	12 (14.63%)	7 (11.29%)
Systolic blood pressure (median) (mm Hg)	129.50 (120.0–143.0)	129.50 (120.0–138.0)
Aneurysm location ^b		
Anterior circulation	77 (96.25%)	59 (98.33%)
Posterior circulation	3 (3.75%)	1 (1.67%)
Aneurysm size (median)	18.00 (13.14–26.0)	17.14 (12.44–24.74)
Aneurysm size classification		
Large (10–15 mm)	30 (36.59%)	24 (38.71%)
Very large or giant (<?>15 mm)	52 (63.41%)	38 (61.29%)

^a Data are No. (%), median (interquartile range), or mean (SD).
^b Of the 6 patients with failed procedures, locations of 2 aneurysms were absent.

Table 3: Sensitivity analyses on the composite primary outcome

Factoring in Subjects Who Missed F/U	6-Mo Occlusion Rate		Unadjusted Result			Adjusted Result ^a		
	Tubridge Group	Control Group	OR	95% CI	P Value	OR	95% CI	P Value
Missing F/U excluded (PPS)	55/73, 75.34%	13/53, 24.53%	9.40	4.14–21.38	<.001	9.31	4.00–21.66	<.001
Missing F/U counted as complete occlusion	64/82, 78.05%	22/62, 35.48%	6.47	3.09–13.52	<.001	6.14	2.92–12.91	<.001
Missing F/U counted as not completely occluded	55/82, 67.07%	13/62, 20.97%	7.68	3.57–16.51	<.001	7.34	3.36–16.05	<.001
Missing F/U counted as complete occlusion for control and incomplete occlusion for Tubridge	55/82, 67.07%	22/62, 35.48%	3.70	1.85–7.42	<.001	3.49	1.73–7.03	<.001

Note:—F/U indicates follow-up; PPS, per-protocol set.
^a Adjusted for each subject's age and aneurysm size (15 versus <?>15 mm).

Diverter (MARCO POLO), FIAT, Flow Diverter Stent for Endovascular Treatment of Unruptured Saccular Wide-necked Intracranial Aneurysms (EVIDENCE), and Complete Occlusion of Coilable Aneurysms (COCO) trials.^{19–22} However, the number of studies that focused on FD treatment of large or giant intracranial aneurysms is low. The purpose of the PARAT trial was to compare the safety and effectiveness outcomes in the treatment of ICA or vertebral artery large or giant aneurysms with the Tubridge FD (with and

without coiling) versus a more conventional approach using Enterprise SAC. The PARAT trial results suggested that in subjects with ICA or vertebral artery large or giant aneurysms, Tubridge FD implantation had a significantly higher 6-month complete occlusion rate compared with conventional stent-assisted coiling. However, there was a trend toward increased risk of stroke with FD implantation.

The treatment of large or giant intracranial aneurysms has evolved significantly during the past few years. Before the emergence of FDs, parent artery occlusion, coiling alone, and stent-assisted coiling were the major treatment modalities for intracranial large or giant aneurysms. According to a meta-

analysis by Turfe et al,²³ parent artery occlusion can result in a complete occlusion rate of 93.0% (95% CI, 86.0% – 97.0%). However, parent artery occlusion may be a viable treatment option only when there is sufficient compensating blood flow. Even when a balloon occlusion test finding is negative, a 4% – 15% complication rate is possible.²³ In addition, there have been concerns about de novo aneurysm occurrence after carotid occlusion.²⁴ In a report by Arambepola et al,²⁵ 4.3% of patients developed de novo aneurysms within a mean of 9 years. Bypass surgery may reduce the incidence of ischemic events, but the procedure may be too complicated, leading to morbidity and mortality rates as high as 7% and 13%, respectively.^{26–28} In this trial, parent artery occlusion was not selected as a feasible control treatment because parent artery sacrifice is considered, throughout most hospitals in China, a last and somewhat futile option for treating intracranial aneurysms. Thus, only those cases posing considerable difficulty or possible failure in parent artery reconstruction would be treated with this method as a salvaging effort. Other than with simple and uncomplicated cases, coiling of large or giant aneurysms without a stent is undertaken far less frequently because of anticipated high recanalization rates.²³

Table 4: Data correlating aneurysm location and size with aneurysm complete occlusion rate

Aneurysm Location	No.	Aneurysm Size (mm)	Proportion of Giant Aneurysms	100% Occlusion Rate at 6-Mo F/U	
				Tubridge	Controls
ICA communicating	19	18.0±7.6 (10.0–33.9)	4/19	76.9%, 10/13	40%, 2/5
ICA ophthalmic	39	15.85.6 (10.0–30.0)	4/39	83.3%, 15/18	33.3%, 6/18
ICA paraclinoid or cavernous	73	21.87.5 (10.0–45)	25/73	75.7%, 28/37	13.8%, 4/29
ICA petrous	3	21.39.3 (15.0–32.0)	1/3	50%, 1/2	0
Vertebral artery	4	14.32.2 (12.0–16.3)	0/4	33.3%, 1/3	100%, 1/1

Note:—F/U indicates follow-up.

Table 5: Secondary outcome evaluation

Secondary Outcome	Enterprise n = 62	Tubridge n = 82	Odds Ratio (95% CI)	P Value
Immediate technique success rate	59/62 (95.16%)	79/82 (96.34%)	1.34 (0.26–6.87)	.726
Death or target vessel-related stroke				
30 days	7/62 (11.29%)	11/82 (13.41%)	0.821 (0.299–2.258)	.703
90 days	8/62 (12.90%)	12/82 (14.63%)	0.864 (0.330–2.263)	.766
1 yr	9/62 (14.52%)	14/82 (17.07%)	0.825 (0.332–2.051)	.678
Mortality rate				
30 days	1/62 (1.61%)	3/82 (3.66%)	0.432 (0.044–4.253)	.634
90 days	1/62 (1.61%)	4/82 (4.88%)	0.320 (0.035–2.934)	.391
1 yr	2/62 (3.23%)	4/82 (4.88%)	0.650 (0.115–3.668)	.699
Hemorrhagic stroke related to target vessel				
30 days	2/62 (3.23%)	5/82 (6.1%)	0.531 (0.096–2.738)	.699
90 days	2/62 (3.23%)	5/82 (6.1%)	0.531 (0.096–2.738)	.699
1 yr	2/62 (3.23%)	5/82 (6.1%)	0.531 (0.096–2.738)	.699
Ischemic stroke related to target vessel				
30 days	4/62 (6.54%)	6/82 (7.32%)	0.874 (0.236–3.239)	.874
90 days	5/62 (8.06%)	6/82 (7.32%)	1.111 (0.323–3.822)	.867
1 yr	6/62 (9.68%)	8/82 (9.76%)	0.991 (0.325–3.109)	.987
Aneurysm rupture				
30 days	1/62 (1.61%)	3/82 (3.66%)	0.357 (0.072–1.783)	.634
90 days	1/62 (1.61%)	3/82 (3.66%)	0.357 (0.072–1.783)	.634
1 yr	1/62 (1.61%)	3/82 (3.66%)	0.357 (0.072–1.783)	.634
General adverse events				
30 days	26/62 (41.94%)	42/82 (51.22%)	0.808 (0.293–2.229)	.269
90 days	26/62 (41.94%)	42/82 (51.22%)	0.808 (0.293–2.229)	.269
1 yr	33/62 (53.23%)	46/82 (56.10%)	0.938 (0.644–1.365)	.88
Rate of intrastent stenosis	2/53 (3.77%)	4/73 (5.48%)	0.676 (0.119–3.837)	1
Rate of intrastent thrombosis	1/59 (1.69%)	6/79 (7.59%)	0.215 (0.025–1.839)	.273

^a Two subjects died during follow-up for reasons unrelated to the target vessels (1 each in Tubridge and control groups).

In many prospective multicenter studies, complete occlusion rates at final follow-up varied from 49% to 93.4%, and 6-month complete occlusion rates ranged from 55.7% to 93.3% (On-line Table 6). In the treatment of large or giant aneurysms with FDs, Becske et al¹⁸ reported a complete occlusion rate of 76.4% at 180 days. Chalouhi et al²⁹ reported a complete occlusion rate of 86% at a median angiographic follow-up of 7 months. On the basis of a meta-analysis of 29 studies, including 1451 patients with 1654 aneurysms, Brinjikji et al³⁰ reported that the complete occlusion rate after FD implantation was 74% (95% CI, 63% – 83%) for large aneurysms and 76% (95% CI, 53% – 90.0%) for giant aneurysms.

Generally, the complete occlusion rate of 75.3% in the Tubridge group at 6 months seems comparable with rates reported in previous publications. The complete occlusion rate of 24.5% for the control group is lower than that reported in previously published articles. In 2 relatively large studies involving large or giant aneurysms, complete occlusion rates of

31% and 41% were reported by Sluzewski et al⁵ and Chalouhi et al,²⁹ respectively. The discrepancies between the latter reports and other published results, as noted above, may be related to variations in patient-selection criteria, study methods, trial design, aneurysm characteristics, and evaluation specifications and procedures. Nevertheless, the effectiveness of Tubridge FD treatment appears to be superior to that of conventional stent-assisted coiling and comparable with reported effectiveness outcomes for other FDs.

Other concerns after FD implantation include technical success rates and overall safety, as well as how FD treatment compares with conventional treatment options. The Enterprise stent is thought to be a safe device for aneurysm treatment, and hemorrhagic complications are assumed to be uncommon.³¹ Although exact causes are not well-understood, delayed aneurysm rupture or intraparenchymal hemorrhage after FD implantation has been recognized.

As shown in On-line Table 6, the hemorrhagic stroke and aneurysm rupture rates were as high as 6.9% and 5.2%, respectively, as reported for multicenter prospective studies. In the International retrospective study of the Pipeline embolization device (IntraPED) study, capturing data from 793 patients among 17 centers, the intraparenchymal hemorrhage rate was 2.4%, while the aneurysm rupture rate was 0.6%.³² In a recent meta-analysis of 3125 treated subjects, the calculated intraparenchymal hemorrhage and aneurysm rupture rates were 2.9% and 1.8%, respectively. These complications are customarily thought to be higher in large or giant aneurysms.³³ Calculated intraparenchymal hemorrhage and aneurysm rupture rates were 5.4% and 7.5%, respectively, in giant aneurysms, and 2.1% and 1.3%, respectively, in small and large aneurysms. The IntraPED study generated similar results, with intraparenchymal hemorrhage and aneurysm rupture rates of 5.8% and 5.8%, respectively, in giant ICA aneurysms; 2.6% and 0.5%, respectively, in large ICA aneurysms; and, 1.9% and 0%, respectively, in small ICA aneurysms. We observed an overall hemorrhagic rate of 6.1% (5/82) in the Tubridge group, of which 2 hemorrhagic occurrences were thought to be the result of procedural injury. Excluding these 2 cases, 2/82 (2.4%) manifested as aneurysm rupture. We did not encounter intraparenchymal hemorrhage, but we observed 1 (1.2%) lateral Sylvian fissure SAH, unrelated to aneurysm rupture. The results of our study are very comparable with those of the above-described studies.

The prospect of ischemic stroke presents additional concerns. A complication rate of 9.68% in the Enterprise control group was consistent with that reported previously by Chalouhi et al.,²⁹ as well as our own single-center experience, in the range of 7.58% – 11.4%.³⁴ Overall, ischemic rates after FD implantation vary among multiple publications, with a rate of 0% – 10.3%. Two recent meta-analyses indicated ischemic rates ranging from 5.5% to 7.5%.^{33,35} However, the occurrence was increased in large or giant aneurysms due to intra-aneurysmal thrombosis or a prolonged procedure time, with a rate ranging from 5.2% to 13.5%.³² Ye et al.³³ indicated an ischemic rate of 9.5% for giant aneurysms. Among Tubridge subjects of the PARAT trial, ischemic complications occurred in 8 (9.76%) subjects. Although the PARAT trial rates are comparable with reported

thromboembolic rates in published studies, a 9% – 10% incidence of ischemic complications should be anticipated when using FDs for large or giant intracranial aneurysms. In the FIAT study, 12/75 subjects (16.0%; 95% CI, 8.9% – 26.7%) treated by flow diversion were dead (n=8) or dependent (n=4) at ≥ 3 months postimplantation.⁸ These results differ from those in our previous experience with FDs. We believe this discrepancy may reflect “real world” early experience with flow diverters for large or giant aneurysms in multiple centers with adjudication by an external imaging Core Laboratory. Similar reports of low procedure-related risks are seen commonly with single-center or retrospective studies, as best exemplified by the recent the Stenting vs. Aggressive Medical Management for Preventing Recurrent Stroke in Intracranial Stenosis (SAMMPRIS) trial.³⁶

To better understand why these complications occurred, we undertook a detailed review of the PARAT trial complications. One possible factor is the role of the learning curve with flow diversion, and physicians’ experience may still play a role in determining the relative safety of FD therapy. The prospect of a learning curve not only encourages a higher skill level in device handling and stent deployment but also leads to more appropriate subject selection. After gaining more experience with a new flow diverter in the treatment of large or giant aneurysms, physicians can avoid many technical failures, make more appropriate FD size selections, better understand the necessity of appropriate postdeployment balloon dilation to avoid poor apposition or overdilation, and better identify the need for simultaneous coil insertion in the treatment of large or giant aneurysms. As an example, the PARAT trial experienced high rates of parent artery occlusions, and these contributed to 3 ischemic events in the Tubridge FD group. The possibility of malapposition may exist when the devices are underdeployed. In such instances, balloon dilation after initial FD deployment may protect against in-stent thrombosis and parent artery occlusion.

The concept of a PARAT trial learning curve was also supported by our multivariate analysis, which showed significantly lower complication rates in the primary investigational site, where the leading physicians already had gained considerable experience with the Tubridge FD in a previous single-center study. We

believe that intensive training of inexperienced physicians on all flow diverters should be advocated before launching a study or introducing use in general neurointerventional practice. We believe that the role that learning curves can play in trials should be increasingly emphasized.

Several recent studies have suggested that antiplatelet regimens may play a role in the occurrence of hemorrhagic or thrombotic complications.^{37,38} In the PARAT trial, investigator management of antiplatelet therapy and subject compliance with prescribed study medications were not investigated.

We have identified some key limitations of the PARAT trial. The most obvious one is the high number of postconsent subject withdrawals (n=41) in both arms, highlighting the delicate balance between physicians and subjects as to what treatment might be better suited for optimal outcomes in the treatment of highly challenging aneurysms. Seven prospective Tubridge subjects voiced their concerns about the potential adverse effects of flow-diverter use, while 15 subjects did not want to undergo stent-assisted coiling because they feared that the long-term durability was inferior to that anticipated with flow diverters. These subjects preferred to withdraw from the trial and wait to be treated with a commercially available flow diverter. The concept of new or apparently complicated technology can be intimidating for some patients.

Among the remaining subjects who departed from this trial prematurely, 1 dropped out after an incomplete subject consent, 6 had their index procedure terminated early after angiography examination indicated safety concerns, 1 was enrolled mistakenly, 1 violated the trial protocol and did not qualify for the index procedure, and 10 were terminated or withdrew from the trial for unknown reasons. Incomplete trial participation by the abovementioned subjects led to a need to expand the sample size and prevent unbalanced subject numbers and age. We compared the baseline data of subjects who were LTF with those from the full analysis set and determined that their removal had not affected the overall distribution of characteristics. Moreover, although subject age was not thought to affect outcomes, we used multivariable logistic regression, adjusted for subject age, to eliminate the effect of this imbalance. Second, the results of this study only

represent the treatment of a specific subtype of aneurysms and a preliminary application of a novel FD. As discussed above, increasing familiarity with the Tubridge device is expected to lower the rate of clinical complications. Finally, although we included treatment of vertebral artery aneurysms in this trial, only 4 subjects with these aneurysms were enrolled. Thus, generalizations about the use of flow diverters in vertebral arteries should be made with caution. Clearly, additional vertebral artery aneurysm studies are warranted.

CONCLUSIONS

This trial showed that there was a significantly higher obliteration rate of aneurysms treated with the Tubridge FD versus conventional Enterprise SAC treatment in a selected group of patients with large or giant intracranial aneurysms. However, this higher obliteration rate comes at a cost of nonsignificant higher complication rates. For inexperienced operators who are either entering flow-diverter clinical trials or preparing for real-world clinical use, the importance of a learning curve should be emphasized.

ACKNOWLEDGMENTS

We thank Professors Ming-hua Li, Qiu-ping Li, and Liang Shao for their contribution in evaluating the 6-month angiographic results; Professors Yin Jian-hua, Chen Yao, Dong Lin, and Dong-Lei Song for data and safety monitoring; and Professors XiaoYan Yan, Chong Ya Dong, and Yao Chen for their statistical analyses. Professor Jia He gave some guidance when drafting and revising this manuscript. We thank them for their constructive guidance.

We are grateful to all principal investigators and delegated physicians who enrolled participants in all participating centers. We also thank the participants and their caretakers and families who consented to participate in the PARAT trial.

REFERENCES

1. International Study of Unruptured Intracranial Aneurysms Investigators. Unruptured intracranial aneurysm: risk of rupture and risks of surgical intervention. *N Engl J Med* 1998;339:1725 – 33 CrossRef Medline
2. Barrow DL, Alleyne C. Natural history of giant intracranial aneurysms and indications for intervention. *Clin Neurosurg* 1995;42: 214 – 44 Medline
3. Darsaut TE, Darsaut NM, Chang SD, et al. Predictors of clinical and angiographic outcome after surgical or endovascular therapy of very large and giant intracranial aneurysms. *Neurosurgery* 2011;68: 903 – 15;

discussion 915 CrossRef Medline

4. Gonzalez NR, Duckwiler G, Jahan R, et al. Challenges in the endovascular treatment of giant intracranial aneurysms. *Neurosurgery* 2008;62(6 Suppl 3):1324 - 35 Medline 5. Sluzewski M, Menovsky T, van Rooij WJ, et al. Coiling of very large or giant cerebral aneurysms: long-term clinical and serial angiographic results. *AJNR Am J Neuroradiol* 2003;24:257 - 62 Medline 6. Chalouhi N, Tjoumakaris S, Gonzalez LF, et al. Coiling of large and giant aneurysms: complications and long-term results of 334 cases. *AJNR Am J Neuroradiol* 2014;35:546 - 52 CrossRef Medline 7. Sadasivan C, Cesar L, Seong J, et al. Treatment of rabbit elastase induced aneurysm models by flow diverters: development of quantifiable indexes of device performance using digital subtraction angiography. *IEEE Trans Med Imaging* 2009;28:1117 - 25 CrossRef Medline 8. Raymond J, Gentric JC, Darsaut TE, et al. Flow diversion in the treatment of aneurysms: a randomized care trial and registry. *J Neurosurg* 2017;127:454 - 62 CrossRef Medline 9. Zhou Y, Yang PF, Fang YB, et al. A novel flow-diverting device (Tubridge) for the treatment of 28 large or giant intracranial aneurysms: a single-center experience. *AJNR Am J Neuroradiol* 2014;35:2326 - 33 CrossRef Medline 10. Kallmes DF, Brinjikji W, Boccardi E, et al. Aneurysm Study of Pipeline in an Observational Registry (ASPIRe). *Interv Neurol* 2016;5: 89 - 99 CrossRef Medline 11. Becske T, Kallmes DF, Saatci I, et al. Pipeline for uncoilable or failed aneurysms: results from a multicenter clinical trial. *Radiology* 2013; 267:858 - 68 CrossRef Medline 12. Byrne JV, Beltechi R, Yarnold JA, et al. Early experience in the treatment of intra-cranial aneurysms by endovascular flow diversion: a multicentre prospective study. *PLoS One* 2010;5. pii: e12492 CrossRef Medline 13. Nelson PK, Lylyk P, Szikora I, et al. The Pipeline embolization device for the intracranial treatment of aneurysms trial. *AJNR Am J Neuroradiol* 2011;32:34 - 40 CrossRef Medline 14. Kan P, Siddiqui AH, Veznedaroglu E, et al. Early postmarket results after treatment of intracranial aneurysms with the Pipeline embolization device: a U.S. multicenter experience. *Neurosurgery* 2012;71: 1080 - 87; discussion 87 - 88 CrossRef Medline 15. Wakhloo AK, Lylyk P, de Vries J, et al; Surpass Study Group. Surpass flow diverter in the treatment of intracranial aneurysms: a prospective multicenter study. *AJNR Am J Neuroradiol* 2015;36:98 107 CrossRef Medline 16. Yu SC, Kwok CK, Cheng PW, et al. Intracranial aneurysms: midterm outcome of Pipeline embolization device—a prospective study in 143 patients with 178 aneurysms. *Radiology* 2012;265:893 - 901 CrossRef Medline 17. Chiu AH, Cheung AK, Wenderoth JD, et al. Long-term follow-up results following elective treatment of unruptured intracranial aneurysms with the Pipeline embolization device. *AJNR Am J Neuroradiol* 2015;36:1728 - 34 CrossRef Medline 18. Becske T, Potts MB, Shapiro M, et al. Pipeline for uncoilable or failed aneurysms: 3-year follow-up results. *J Neurosurg* 2017;127:81 - 88 CrossRef Medline 19. ClinicalTrials.gov. Multicenter Randomized Trial on Selective Endovascular Aneurysm Occlusion With Coils Versus Parent Vessel Reconstruction Using the SILK Flow Diverter (MARCO POLO Post Market Clinical Investigation). <https://clinicaltrials.gov/ct2/show/NCT01084681>. Accessed March 17, 2014 20. ClinicalTrials.gov. A Randomized Trial Comparing Flow Diversion and Best-standard Treatment: the FIAT Trial. <https://www.clinicaltrials.gov/ct2/show/NCT01349582>. Accessed January 20, 2016 21. ClinicalTrials.gov. Multicenter Randomized Study for Medico economic Evaluation of Embolization With Flow Diverter Stent in the Endovascular Treatment of Unruptured Saccular Wide-necked Intracranial Aneurysms. <https://www.clinicaltrials.gov/ct2/show/NCT01811134>. Accessed January 20, 2016

22. ClinicalTrials.gov. Complete Occlusion of Coilable Aneurysms. <https://www.clinicaltrials.gov/ct2/show/NCT00777907>. Accessed January 20, 2016 23. Turfe ZA, Brinjikji W, Murad MH, et al. Endovascular coiling versus parent artery occlusion for treatment of cavernous carotid aneurysms: a meta-analysis. *J Neurointerv Surg* 2015;7:250 - 55 CrossRef Medline 24. Zhou Y, Yang P, Zhang Y, et al. Posterior cerebral artery-posterior communicating artery (PCA-PComA) aneurysms: report of five cases and literature review. *Neurol India* 2012;60:228 - 30 CrossRef Medline 25. Arambepola PK, McEvoy SD, Bulsara KR. De novo aneurysm formation after carotid artery occlusion for cerebral aneurysms. *Skull Base* 2010;20:405 - 08 CrossRef Medline 26. Jafar JJ, Russell SM, Woo HH. Treatment of giant intracranial aneurysms with saphenous vein extracranial-to-intracranial bypass grafting: indications, operative technique, and results in 29 patients. *Neurosurgery* 2002;51:138 - 44; discussion 44 - 46 CrossRef Medline 27. Spetzler RF, Schuster H, Roski RA. Elective extracranial intracranial arterial bypass in the treatment of inoperable giant aneurysms of the internal carotid artery. *J Neurosurg* 1980;53:22 - 27 CrossRef Medline 28. Morgan MK, Ferch RD, Little NS, et al. Bypass to the intracranial internal carotid artery. *J Clin Neurosci* 2002;9:418 - 24 CrossRef Medline 29. Chalouhi N, Tjoumakaris S, Starke RM, et al. Comparison of flow diversion and coiling in large unruptured intracranial saccular aneurysms. *Stroke* 2013;44:2150 - 54 CrossRef Medline 30. Brinjikji W, Murad MH, Lanzino G, et al. Endovascular treatment of intracranial aneurysms with flow diverters: a meta-analysis. *Stroke* 2013;44:442 - 47 CrossRef Medline 31. Fargen KM, Hoh BL, Welch BG, et al. Long-term results of Enterprise stent-assisted coiling of cerebral aneurysms. *Neurosurgery* 2012;71:239 - 44; discussion 244 CrossRef Medline 32. Kallmes DF, Hanel R, Lopes D, et al. International retrospective study of the Pipeline embolization device: a multicenter aneurysm treatment study. *AJNR Am J Neuroradiol* 2015;36:108 - 15 CrossRef Medline 33. Ye G, Zhang M, Deng L, et al. Meta-analysis of the efficiency and prognosis of intracranial aneurysm treated with flow diverter devices. *J Mol Neurosci* 2016;59:158 - 67 CrossRef Medline 34. Zhang Y, Zhou Y, Yang P, et al. Comparison of the flow diverter and stent-assisted coiling in large and giant aneurysms: safety and efficacy based on a propensity score-matched analysis. *Eur Radiol* 2016; 26:2369 - 77 CrossRef Medline 35. Zhou G, Su M, Yin YL, et al. Complications associated with the use of flow-diverting devices for cerebral aneurysms: a systematic review and meta-analysis. *Neurosurg Focus* 2017;42:E17 CrossRef Medline 36. Derdeyn CP, Chimowitz MI, Lynn MJ, et al; Stenting and Aggressive Medical Management for Preventing Recurrent Stroke in Intracranial Stenosis Trial Investigators. Aggressive medical treatment with or without stenting in high-risk patients with intracranial artery stenosis (SAMMPRIS): the final results of a randomised trial. *Lancet* 2014;383:333 - 41 CrossRef Medline 37. Delgado Almandoz JE, Crandall BM, Scholz JM, et al. Last recorded P2Y12 reaction units value is strongly associated with thromboembolic and hemorrhagic complications occurring up to 6 months after treatment in patients with cerebral aneurysms treated with the Pipeline embolization device. *AJNR Am J Neuroradiol* 2014;35: 128 - 35 CrossRef Medline 38. Skukalek SL, Winkler AM, Kang J, et al. Effect of antiplatelet therapy and platelet function testing on hemorrhagic and thrombotic complications in patients with cerebral aneurysms treated with the Pipeline embolization device: a review and meta-analysis. *J Neurointerv Surg* 2016;8:58 - 65 CrossRef Medline 39. McAuliffe W, Wycoco V, Rice H, et al. Immediate and midterm results following treatment of unruptured intracranial aneurysms with the pipeline embolization device. *AJNR Am J Neuroradiol* 2012; 33:164 - 70 CrossRef Medline 10 Liu

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Health Economic Evaluation on Domestic and Imported Devices of Radiofrequency Ablation in Persistent Atrial Fibrillation

Li Xiang et al

Pecking University Clinical Research Institute, Beijing, 100871

Abstract Objective: To analyze cost utility between Columbus™ and Carto™ 3 from the payer's perspective. Methods: This study developed a one-year decision model to compare the health outcome and cost. The cost data of only medical devices cost was collected from payer's perspective. The sensitivity analysis was conducted. Results: The health outcome were 0.7647 QALY and 0.7681 QALY separately for domestic and imported devices. The cost per surgical operation were 42566.34 RMB and 66332 RMB

for domestic and imported devices. Compared to domestic devices, 6989900 RMB would be cost for one more QALY gained, which was much more than the willingness-to-pay of 161940 RMB. Conclusion: This research showed that domestic devices had better cost utility than imported devices, which is proposed to be preferable considered compared to imported devices in decision making. Key Words Radiofrequency Ablation; Persistent Atrial Fibrillation; Cost Utility Analysis

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- [1] GO A S. Prevalence of diagnosed atrial fibrillation in adults: national implications for rhythm management and stroke prevention: the AnTicoagulation and Risk Factors in Atrial Fibrillation (ATRIA) Study[J]. JAMA, 2001,285(18):2370–2375.
- [2] 黄从新. 中国经导管消融治疗心房颤动注册研究 10 年回顾 [J]. 中华心律失常学杂志, 2009,13(4):296–302.
- [3] WILBER D J. Comparison of antiarrhythmic drug therapy and radiofrequency catheter ablation in patients with paroxysmal atrial fibrillation: a randomized controlled trial[J]. JAMA, 2010,303(4): 333–340.
- [4] 上海微创电生理医疗科技有限公司, 临床试验总结报告, 2015:501–503.
- [5] POKUSHALOV E. Ablation of paroxysmal and persistent atrial fibrillation: 1-year follow-up through continuous subcutaneous monitoring [J]. J Cardiovasc Electrophysiol,2011,22(4): 369–375.
- [6] 鲁志兵. 国产 ColumbusTM 三维标测系统在心房颤动射频消融中的初步应用 [J]. 中国心脏起搏与心电生理杂志, 2014(6):532–534.
- [7] 金志清. 国产心脏三维电解剖标测系统引导下射频导管消融持续性心房颤动 [J]. 中华心律失常学杂志, 2015, 19(4): 275–278.
- [8] 常栋. 新型三维心脏电生理标测系统 ColumbusTM 在心房颤动导管消融中的应用 [J]. 中华生物医学工程杂志, 2015, 21(5): 441–444.
- [9] KETTERING K, YIM D H, GRAMLEY F. Catheter ablation of persistent atrial fibrillation: circumferential pulmonary vein ablation: beneficial effect of an additional linear lesion at the roof of the left atrium on the long-term outcome[J]. Herzschrittmacherther elektrophysiol, 2017, 28 (3):328–334.
- [10] 上海微创电生理医疗科技有限公司, 微创电生理三维心脏电生理标测系统及冷盐水灌注射频消融导管创新器械上市效应综合评价, 2017.
- [11] KIMURA T. A cost-utility analysis for catheter ablation of atrial fibrillation in combination with warfarin and dabigatran based on the CHADS2 score in Japan[J]. J Cardiol, 2017, 69(1): 89–97.
- [12] JAIS P. Catheter ablation versus antiarrhythmic drugs for atrial fibrillation: the A4 study[J]. Circulation, 2008, 118(24): 2498–2505.
- [13] REYNOLDS M R. Health state utilities improve with catheter ablation for AF more than with drug therapy: Analysis from a randomized trial. 2007.
- [14] 黄从新, 中国房颤导管消融现状 [J]. 医学与哲学, 2016, 37(11): 10–13.

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文 | Diana York

What 'being happy' means to me:

Getting up in the morning, and not dreading going to work, but actually looking forward to the day ahead.

Working for a company like MicroPort whose culture is the reason I look forward to going to work each day.

It is a great feeling to work with people who care about their job and their coworkers. People that are proud of what they do and know their work is valued, recognized and appreciated!

I work with teams in DC at MPO we are like family. Everyone treats each other with respect. Team members are willing to work wherever, whenever to get the job done.

Everyone truly cares about those they work with and it is evident day in and day out.

While salaries, bonuses, vacation, are all great things, knowing that our work contributes to making life better for someone, is most rewarding! How many people can say their company/job improves someone's life!

Great people, great company, these are the things that make me happy in my job! 🇺🇸

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132、Many of the great changes in our lives can be traced to small things, a chance

acquaintance, an accident, or some little happening. A time comes to every man when he must do a thing or miss a great benefit. If the man does it, all is well. If not, it isn't likely that he will have the chance again. You can call that opportunity if you wish, but it is foresight that leads a man to take advantage of the condition of things. Foresight is a most valuable thing to have. —Theodore Roosevelt

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Vision

Bring About a People-Centric SuperCluster of Emerging Medical Technologies

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Mission

Provide the Best Yet Accessible Ways & Means to Prolong and Reshape Lives

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Faith

Help Live Beyond 115 Years of Very Healthy Life for Billions of People on Earth

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Branding Ideology

MicroPort® - The Patient Always Comes First

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Values

Quality Integrity Accountability Efficiency Innovation Aspiration Dedication Collaboration

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Going to the Boundary Waters is the perfect trip for those looking to embrace wilderness and escape society for a few days. Located on the border (hence the name) of the Canada and the US state, Minnesota, the Boundary Waters are a unique landscape in the world. Its passageway of lakes and ponds were carved by glaciers millions of years ago, creating a maze of waterways. Paired with Quetico (the Canadian part of the Waters), the area stretches over 9,000 square kilometers.

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Because the Boundary Waters are part of the federally protected Superior National Forest, there is next to no human activity in the area, save for a limited number of reserved spaces for campers looking for a pristine wilderness experience.

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What should one expect to do in the Boundary Waters? Before even entering the area, campers must prepare for their trek. The remoteness of the area means that prep work is extremely important, and forgetting even one tool could severely hinder one's stay. A basic check list includes: camping permit, canoes, pots and pans, utensils, bowls, food, lighter, camping packs, tent, map, compass, water purifiers, clothes, flashlight, sleeping bag, sunscreen, and mosquito repellent (lots of it). Without these items, it is impossible to camp in the boundary waters. These are necessities. Some other items I personally recommend bringing are fire-starters (cotton balls soaked in wax melted in a microwave), a portable hammock, a book, playing cards, and a pair of extra comfy socks that should be kept dry at all times!

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
Once the above items have been rounded up, you must plan your route through the maze of lakes. Because there is no cell coverage, the map and compass are your only tools for navigation. Camp sites and “portages” are found on the maps to help guide and orient campers. These portage paths refer to the spots on lakes where campers de-canoe and haul their supplies to the next lake on foot. Most of these are relatively short treks, between 100 and 300 meters, but others can be up to 3 kilometers long. Imagine carrying 50 kilograms of supplies, including a canoe on your head, for that long!

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Though the Boundary Waters may not be suitable for the beginning camper, those who are comfortable in the wilderness will find serenity and beauty seldom seen anywhere else.

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My favorite part of each day is setting up camp, cooking food, and relaxing with friends before another day on the lakes. Most days we paddle between 10 and 15 kilometers, and are exhausted and hungry by the time we finish setting up our tent. A warm meal, followed by hilarious conversations, fun games we make up on the spot, and cool evening swims are unforgettable. As is the view of the Milky Way galaxy. One night I counted over 50 meteors fall on a background of thousands of bright stars, some exploding before our eyes as they crashed down to earth hundreds of kilometers away.

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Camping in the Boundary Waters is an eye-opening experience, as one removes themselves from society and reflects on life between each stroke of the paddle. There is much wisdom to be found in the woods and on the lakes. Just be sure to bring plenty of mosquito repellent. They're called the Minnesota State Bird for a reason!

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